

Anti-TLR2 antibody [mT2.7] (STJ16100485)

STJ16100485

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR2 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunoprecipitation and Immunohistochemistry research applications.
Applications	IHC-F/FC/ELISA/IP/IHC-P
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	mT2.7
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7097
Gene Symbol	TLR2
Uniprot ID	TLR2_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	