

Anti-TLR4 antibody [HTA125] (STJ16100496)

STJ16100496

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR4 is suitable for use in Flow Cytometry, ELISA, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	FC/ELISA/IF/IP/WB
Host/Source	Mouse
Reactivity	Human/Dog/Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HTA125
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7099
Gene Symbol	TLR4
Uniprot ID	TLR4_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	