

Anti-CD32 antibody [AT10] {Biotin} (STJ16100503)

STJ16100503

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD32 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human/Pig/Rhesus/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AT10
Concentration	0.1 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2212
Gene Symbol	FCGR2A
Uniprot ID	FCG2A_HUMAN
Immunogen	K562 cell line.
Immunogen Region	
Specificity	
Immunogen Sequence	