

Anti-CD62L antibody [FMC46] {FITC} (STJ16100506)

STJ16100506

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD62L is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human/Rhesus/Bovine/Cynomolgus Monkey/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FMC46
Concentration	0.1 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	PHA stimulated lymphoblasts
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081