

Anti-CD102 antibody [B-T1] (STJ16100514)

STJ16100514

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD102 is suitable for use in Western Blot, Immunofluorescence, Flow Cytometry and Immunoprecipitation research applications.
Applications	WB/IF/FC/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B-T1
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	ICAM-2 transfected CHO cells.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	