

## Anti-CD104 antibody [450-9D] (STJ16100525)

STJ16100525

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD104 is suitable for use in Western Blot, Immunoprecipitation, Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	WB/IP/FC/IHC-F/IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Ferret/Mink/Mustelid

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	450-9D
<b>Concentration</b>	1.0 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Purified alpha6 beta4 integrin from A431 cells.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	

