

Anti-CD152 antibody [BNI3] (STJ16100529)

STJ16100529

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD152 is suitable for use in ELISA, Immunofluorescence, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	ELISA/IF/FC/IHC-F/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BNI3
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	1493
Gene Symbol	CTLA4
Uniprot ID	CTLA4_HUMAN
Immunogen	Human CTLA-4/human IgG heavy chain fusion protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	