

Anti-CD89 antibody [MIP8a] (STJ16100532)

STJ16100532

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD89 is suitable for use in Immunoprecipitation, Flow Cytometry and ELISA research applications.
Applications	IP/FC/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MIP8a
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2204
Gene Symbol	FCAR
Uniprot ID	FCAR_HUMAN
Immunogen	Recombinant soluble human Fc alpha R.
Immunogen Region	
Specificity	
Immunogen	
Sequence	