

Anti-CD91 antibody [A2Mr alpha-2] (STJ16100538)

STJ16100538

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD91 is suitable for use in Western Blot, Immunoprecipitation, Flow Cytometry, Immunohistochemistry and ELISA research applications.
Applications	WB/IP/FC/IHC-F/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A2Mr alpha-2
Concentration	1.0mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified alpha2 macroglobulin receptor
Immunogen Region	
Specificity	
Immunogen Sequence	