

Anti-CD239 antibody [BRIC221] (STJ16100542)

STJ16100542

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD239 is suitable for use in Western Blot, ELISA, Immunofluorescence, Flow Cytometry and Immunohistochemistry research applications.
Applications	WB/ELISA/IF/FC/IHC-F
Host/Source	Mouse
Reactivity	Human/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BRIC221
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	TRIS buffered glycine
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human erythrocytes.
Immunogen Region	
Specificity	
Immunogen Sequence	