

Anti-CD173 antibody [BRIC231] (STJ16100544)

STJ16100544

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD173 is suitable for use in Immunofluorescence and Flow Cytometry research applications.
Applications	IF/FC
Host/Source	Mouse
Reactivity	Human/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BRIC231
Concentration	1.0mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	TRIS buffered glycine
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human erythroleukemic cell line (HEL) established from a 30 year old patient with relapsed erythroleukemia following treatment for Hodgkin lymphoma.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	