

Anti-CD175 antibody [BRIC111] (STJ16100545)

STJ16100545

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD175 is suitable for use in ELISA and Immunohistochemistry research applications.
Applications	ELISA/IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BRIC111
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	TRIS buffered glycine
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human Tn red blood cells.
Immunogen Region	
Specificity	
Immunogen Sequence	