

**Anti-CD88 antibody [P12/1] (STJ16100547)**

STJ16100547

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD88 is suitable for use in Western Blot, Immunohistochemistry, Flow Cytometry and Immunohistochemistry research applications.
<b>Applications</b>	WB/IHC-P/FC/IHC-F
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Rhesus

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	P12/1
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgG2a
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	C5aR-peptide: Met1-Asn31
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	