

Anti-CD146 antibody [OJ79c] (STJ16100552)

STJ16100552

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD146 is suitable for use in Flow Cytometry, Immunohistochemistry, ELISA and Immunofluorescence research applications.
Applications	FC/IHC-F/ELISA/IF
Host/Source	Mouse
Reactivity	Human/Pig/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OJ79c
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant human MUC18 (D1-D5) Fc protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	