

Anti-CD146 antibody [OJ79c] {FITC} (STJ16100553)
STJ16100553**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD146 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human/Pig/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OJ79c
Concentration	0.1 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant human MUC18 (D1-D5) Fc protein.
Immunogen Region	
Specificity	
Immunogen Sequence	