

**Anti-CD195 antibody [HEK/1/85a] (STJ16100555)**  
STJ16100555**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-CD195 is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	HEK/1/85a
<b>Concentration</b>	1.0 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgG2a
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	CHO cells transfected with human CCR5
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	