

Anti-CD195 antibody [HEK/1/85a] {FITC} (STJ16100556)
STJ16100556**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD195 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HEK/1/85a
Concentration	0.1 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	CHO cells transfected with human CCR5
Immunogen Region	
Specificity	
Immunogen Sequence	