

Anti-CD79 antibody [AT107-2] (STJ16100558)

STJ16100558

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD79 is suitable for use in Western Blot, Immunoprecipitation, Flow Cytometry and ELISA research applications.
Applications	WB/IP/FC/ELISA
Host/Source	Rat
Reactivity	Human/Pig/Rat/Tasmanian Devil/Dog/Cat/Human/Koala

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AT107-2
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	TBS with 0, 09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Peptide containing 20 amino acid residues from mouse CD79b conjugated to Keyhole Limpet Hemocyanin (KLH).
Immunogen Region	
Specificity	
Immunogen	
Sequence	