

Anti-CD81 antibody [1D6] (STJ16100575)

STJ16100575

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD81 is suitable for use in Western Blot, Immunohistochemistry, Immunoprecipitation, Flow Cytometry and Immunofluorescence research applications.
Applications	WB/IHC-P/IP/FC/IF
Host/Source	Mouse
Reactivity	Human/Sheep/Chimpanzee/Goat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1D6
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	975
Gene Symbol	CD81
Uniprot ID	CD81_HUMAN
Immunogen	OCI-LY8 cells aggregated by 5A6 (another CD81 antibody)
Immunogen Region	
Specificity	
Immunogen Sequence	