

Anti-kappa light chain antibody [2B7] (STJ16100578)

STJ16100578

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-kappa light chain is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2B7
Concentration	10 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F 1:30 WB
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	kappa Bence Jones protein
Immunogen Region	
Specificity	
Immunogen Sequence	