

Anti-lambda light chain antibody [48] (STJ16100579)

STJ16100579

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-lambda light chain is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-P/IHC-F/WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	48
Concentration	15 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:50 IHC-F 1:50 WB ELISA
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Bence Jones lambda protein
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081