

**Anti-Immunoglobulin IgA antibody [MH14/1] (STJ16100584)**

STJ16100584

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Immunoglobulin IgA is suitable for use in Immunohistochemistry, ELISA and Immunofluorescence research applications.
<b>Applications</b>	IHC-F/IHC-P/ELISA/IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	MH14/1
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC-F IHC-P 1:5-1:10 ELISA IF
<b>Formulation</b>	Contains 0.09% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Isolated human IgA
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	