

Anti-Immunoglobulin IgE antibody [MH25/1] (STJ16100585)

STJ16100585

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Immunoglobulin IgE is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-F/IHC-P/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MH25/1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P 1:5-1:10 IF
Formulation	PBS with 0.5% BSA and 0.09% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Isolated human immunoglobulins
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	