

Anti-Interleukin 6 antibody [8F4] (STJ16100598)

STJ16100598

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Interleukin 6 is suitable for use in FUNC, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	FUNC/ELISA/IP/WB
Host/Source	Mouse
Reactivity	Human/Cynomolgus Monkey/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	8F4
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3569
Gene Symbol	IL6
Uniprot ID	IL6_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	