

Anti-Insulin antibody [IN-05] (STJ16100599)

STJ16100599

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Insulin is suitable for use in ELISA, Immunohistochemistry and Immunocytochemistry research applications.
Applications	ELISA/IHC-P/ICC
Host/Source	Mouse
Reactivity	Pig/Cow/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IN-05
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Porcine insulin.
Immunogen Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081