

Anti-IFN-gamma antibody [F1] (STJ16100652)

STJ16100652

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IFN-gamma is suitable for use in ELISA and Immunoprecipitation research applications.
Applications	ELISA/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	F1
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3458
Gene Symbol	IFNG
Uniprot ID	IFNG_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	