

Anti-IL-6 antibody [MP5-20F3] (STJ16100654)

STJ16100654

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IL-6 is suitable for use in FUNC, ELISA and Immunoprecipitation research applications.
Applications	FUNC/ELISA/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MP5-20F3
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3569
Gene Symbol	IL6
Uniprot ID	IL6_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	