

Anti-Interleukin-1 receptor/P80 antibody [D14e3] (STJ16100658)
STJ16100658**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-1 receptor/P80 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	D14e3
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence