

## Anti-IFN-gamma antibody [F14] (STJ16100661)

STJ16100661

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-IFN-gamma is suitable for use in Immunohistochemistry, ELISA and Immunoprecipitation research applications.
<b>Applications</b>	IHC-F/ELISA/IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	F14
<b>Concentration</b>	100 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3458</a>
<b>Gene Symbol</b>	<a href="#">IFNG</a>
<b>Uniprot ID</b>	<a href="#">IFNG_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	