

Anti-CD200R antibody [OX108] (STJ16100671)

STJ16100671

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD200R is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX108
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	131450
Gene Symbol	CD200R1
Uniprot ID	MO2R1_HUMAN
Immunogen	Fusion protein hCD200RCD4d3+4.
Immunogen Region	
Specificity	
Immunogen Sequence	