

Anti-CD90 antibody [F15-42-1] (STJ16100673)

STJ16100673

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD90 is suitable for use in Immunofluorescence, Immunoprecipitation, Flow Cytometry and Immunohistochemistry research applications.
Applications	IF/IP/FC/IHC-F
Host/Source	Mouse
Reactivity	Human/Cynomolgus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	F15-42-1
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified human brain Thy-1.
Immunogen Region	
Specificity	
Immunogen Sequence	