

Anti-Liver Fatty Acid Binding Protein antibody [L2B10] (STJ16100686)

STJ16100686

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Liver Fatty Acid Binding Protein is suitable for use in Immunohistochemistry, ELISA, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/ELISA/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Baboon/Dog/Rat/Swine

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	L2B10
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	