

Anti-Lyn antibody [LYN-01] (STJ16100692)

STJ16100692

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lyn is suitable for use in Immunoprecipitation, Western Blot and Immunocytochemistry research applications.
Applications	IP/WB/ICC
Host/Source	Mouse
Reactivity	Rat/Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LYN-01
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4067
Gene Symbol	LYN
Uniprot ID	LYN_HUMAN
Immunogen	Bacterially expressed recombinant fragment of human Lyn (aa 8-238).
Immunogen Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081