

**Anti-NHERF antibody [EBP-10] (STJ16100695)**

STJ16100695

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-NHERF is suitable for use in Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
<b>Applications</b>	IP/WB/IHC-P
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	EBP-10
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
<b>Isotype</b>	IgG1k
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	9368
<b>Gene Symbol</b>	NHERF1
<b>Uniprot ID</b>	NHRF1_HUMAN
<b>Immunogen</b>	Bacterially produced recombinant full-length human NHERF1.
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	