

Anti-Matrix Metalloproteinase-2 antibody [4D3] (STJ16100697)

STJ16100697

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Matrix Metalloproteinase-2 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4D3
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4313
Gene Symbol	MMP2
Uniprot ID	MMP2_HUMAN
Immunogen	human MMP-2 peptide: TSLGLPPDVQRVD
Immunogen Region	
Specificity	
Immunogen Sequence	