

Anti-MUC1/Polymorphic epithelial mucin antibody [2G7] (STJ16100705)

STJ16100705

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MUC1/Polymorphic epithelial mucin is suitable for use in ELISA, Immunohistochemistry, Immunofluorescence and Flow Cytometry research applications.
Applications	ELISA/IHC-F/IHC-P/IF/FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2G7
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	4582
Gene Symbol	MUC1
Uniprot ID	MUC1_HUMAN
Immunogen	VNTR 60mer glycosylated synthetic peptide BSA conjugate
Immunogen Region	
Specificity	
Immunogen Sequence	

