

Anti-MUC5AC/M1 antibody [2-11M1] (STJ16100707)

STJ16100707

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MUC5AC/M1 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Cat/Cow/Human/Monkey/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2-11M1
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT). Immunoblotting (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; staining of formalin-fixed tissues requires 0, 1% pronase digestion (0, 1% CaCl ₂ , 0, 1N NaOH) for
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	4586
Gene Symbol	MUC5AC
Uniprot ID	MUC5A_HUMAN
Immunogen	A BALB/c mouse was immunized with mucin isolated from an ovarian cyst fluid (pure endocervical type according to the Fenoglio's classification) from an ALab patient. Fusion partner: Sp2/0.
Immunogen Region	
Specificity	
Immunogen Sequence	

