

Anti-MMP-3 antibody [1B4] (STJ16100712)

STJ16100712

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MMP-3 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P/IHC-F/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1B4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4314
Gene Symbol	MMP3
Uniprot ID	MMP3_HUMAN
Immunogen	human MMP-3 peptide: KSLRKLEPELH
Immunogen Region	
Specificity	
Immunogen Sequence	