

Anti-CD45 antibody [ANK74] (STJ16100717)

STJ16100717

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD45 is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-F/FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ANK74
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	5788
Gene Symbol	PTPRC
Uniprot ID	PTPRC_HUMAN
Immunogen	IL-2-activated cultured NK cells of Wag rats
Immunogen Region	
Specificity	
Immunogen Sequence	