

Anti-T-cell Receptor antibody [R73-2b] (STJ16100718)

STJ16100718

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-T-cell Receptor is suitable for use in Immunohistochemistry and Flow Cytometry research applications. |
| Applications | IHC-F/FC |
| Host/Source | Mouse |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | R73-2b |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | Culture supernatant with 0.05% sodium azide |
| Isotype | IgG2b |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| |
|-------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen |
| Sequence |