

Anti-Hepatitis B virus antibody [HB3] (STJ16100722)

STJ16100722

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Hepatitis B virus is suitable for use in Immunocytochemistry and ELISA research applications.
Applications	ICC/ELISA
Host/Source	Mouse
Reactivity	HBV

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HB3
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified HBsAg from human plasma.
Immunogen Region	
Specificity	
Immunogen Sequence	