

Anti-Laminin antibody [MEC-5] (STJ16100725)

STJ16100725

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Laminin is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Rat
Reactivity	Rat/Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEC-5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	auto-antibody
Immunogen Region	
Specificity	
Immunogen Sequence	