

Anti-Tumor Antigens of Epithelial Origin antibody [CC52] (STJ16100726)

STJ16100726

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tumor antigens of Epithelial Origin is suitable for use in Immunohistochemistry, Flow Cytometry, Western Blot, Immunoprecipitation and ELISA research applications.
Applications	IHC-P/IHC-F/FC/WB/IP/ELISA
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CC52
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	rat colon adenocarcinoma cells (CC531)
Immunogen Region	
Specificity	
Immunogen Sequence	