

Anti-Influenza A Nucleoprotein antibody [AA5H] (STJ16100728)

STJ16100728

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Influenza A Nucleoprotein is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	WB/IHC-P/IF
Host/Source	Mouse
Reactivity	Virus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AA5H
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Influenza A/Puerto Rico/8/34 (H1N1) and A/Bangkok/1/79 (H3N2) viruses.
Immunogen Region	
Specificity	
Immunogen Sequence	