

Anti-Melanoma Discr. Antigen antibody [PAL-M2] (STJ16100731)

STJ16100731

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Melanoma Discr. antigen is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PAL-M2
Concentration	10 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F 1:20
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	human melanoma metastases tissue
Immunogen Region	
Specificity	
Immunogen Sequence	