

Anti-B-cells antibody [MB2] (STJ16100750)

STJ16100750

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-B-cells is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MB2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P 1:5-1:40 WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Isolated from human lymphocytes
Immunogen Region	
Specificity	
Immunogen Sequence	