

Anti-MCP-1 antibody [MNA.1] (STJ16100753)

STJ16100753

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MCP-1 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/FC/ELISA/IP/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MNA.1
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1 % BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	6347
Gene Symbol	CCL2
Uniprot ID	CCL2_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	