

Anti-IP-10 antibody [6D4] (STJ16100756)

STJ16100756

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IP-10 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunofluorescence and Western Blot research applications.
Applications	IHC-F/FC/ELISA/IF/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6D4
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	3627
Gene Symbol	CXCL10
Uniprot ID	CXL10_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	