

Anti-MGL antibody [ER-MP23] (STJ16100767)

STJ16100767

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MGL is suitable for use in Immunohistochemistry, Flow Cytometry, FUNC and Immunohistochemistry research applications.
Applications	IHC-F/FC/FUNC/IHC-P
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ER-MP23
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	11343
Gene Symbol	MGLL
Uniprot ID	MGLL_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	