

Anti-BPI antibody [4H5] (STJ16100769)

STJ16100769

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-BPI is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4H5
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	671
Gene Symbol	BPI
Uniprot ID	BPI_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081